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## Product Description

Medicvent's double mask is a unique reusable anaesthetic face mask that combines two masks in one: a soft inner silicone mask that administers medical gases to the patient; and an outer polysulfone mask that allows it to capture leaked gases from the mask for further evacuation. The leaked gases from the mask can be aspirated by a central fan unit and evacuated out of the room, reducing the concentration of waste gases up to 92% compared to traditional face masks.

The double mask is also available with a neckband, which allows the mask to remain in the correct position for gas evacuation even when the patient is not deliberately holding it.

The double mask is used in combination with Medicvent's coupling house or demand coupling house (accessories), that connects the double mask with a standard anaesthesia breathing tube or analgesic demand valve; and a Ø28 mm scavenging tube.

## Intended purpose

The double mask with accessories is intended to administer medical gases to patients for sedation and short-term pain relief during medical procedures and simultaneously enable evacuation of leaked gases from the mask.

## Patient target group and users

To be used on adults and children older than 1.

Only to be used by healthcare professionals or patients under the supervision of healthcare professionals.

## Indications and clinical benefits

The double mask provides a non-invasive means to administer medical gases for short-term pain relief or sedation into the patient's mouth and nose and decreases exposure of patients and users to excess anaesthetic gases.

It is indicated when inhaled anaesthesia/analgesia is required during medical procedures performed in a hospital and other healthcare environments.

The double mask with coupling house is used in combination with a standard anaesthetic breathing system.

The double mask with demand coupling house is used in combination with a BPR (Carnét) or O\_two analgesic demand valve.

## Contraindications and adverse events

There are no known contraindications or adverse events.

**Any serious incident that has occurred in relation to the devices must be reported to Medicvent AB and the competent authority of the country in which the user and/or the patient is established.**

## Size selection

The double mask is available in the standard sizes 1-3 for children and 4-5 for adults. Choose a size that ensures a complete seal whilst exerting minimum pressure.

## Instructions for Use

- 1 Ensure that the double mask and accessories have been cleaned and disinfected prior to initial use and between patients. See the Reprocessing Instructions.
- 2 Inspect the devices and check the expiry date. The date shall be less than 5 years from the manufacturing date (page 3, figure 1). Discard damaged or deformed material and expired products.  
**WARNING!** Do not use if the devices are cracked or damaged  
**WARNING!** Do not use after expiry date
- 3 Check that the O-ring in the coupling house or demand coupling house is moistened with Fomblin lubricant or a similar lubricant.
- 4 Ensure the correct coupling house or demand coupling house is used.
- 5 As shown in figure 2, page 3, connect a Ø28 mm scavenging tube to the side entry port of the coupling house or demand coupling house, and a Ø22 mm breathing tube or a demand valve to the top entry port.
- 6 Connect the double mask to the bottom entry port of the coupling or demand coupling house.
- 7 Ensure that all components are properly attached, and do not detach.
- 8 In case of a double mask with neckband: Ensure that the neckband is securely attached around the patient's neck.
- 9 Initiate the gas delivery and suction according to manufacturer instructions.
- 10 After completing the procedure, disconnect the double mask and other tubing from the coupling house and demand coupling house by gently twisting and pulling them out.
- 11 Clean and disinfect according to reprocessing instructions.

## Compatibility with other devices

- The coupling house is designed to connect a standard Ø22mm breathing tube and a Ø28 mm scavenging tube.
- The demand coupling houses are designed to connect a demand valve Carnét or O\_two and a Ø28 mm scavenging tube.
- For optimal gas evacuation, use a suction flow between 27 m<sup>3</sup>/h and 35 m<sup>3</sup>/h.

**WARNING!** Do not connect to systems intended for plume evacuation, AGSS or vacuum.

**WARNING!** MRI compatibility has not been tested.

## Disposal

Disposal must be done according to hospital guidelines.

## General Recommendations

- As the double mask enables the evacuation of any leaking gas when connected to a local or central evacuation system, there is no need to hold the mask tightly against the face or to turn off/decrease the gas flow when the double mask is removed or pulled aside.

- During surgery, the double mask is placed on the patient's face with a light pressure with the fingertips (thumb and index finger) to keep the mask in place. The palm of the hand and other fingers may be used to reposition the patient's chin to ensure that the airways are free.

## Reprocessing instructions

<b>Manufacturer:</b> Medicvent AB <b>Processing Methods:</b> Cleaning, Disinfection, Drying		  
<b>Processing instructions for:</b> <ul style="list-style-type: none"> <li>• 210001-210005 Double Mask, Size 1 to 5</li> <li>• 210013-210015 Double Mask with neckband, Size 3 to 5</li> <li>• 210006            Coupling house</li> <li>• 130006            Demand O_two Coupling House</li> <li>• 130007            Demand Carnét Coupling House</li> </ul>		
<b>WARNINGS</b>	<ul style="list-style-type: none"> <li>• Do not use if the devices are cracked or damaged</li> <li>• Do not use after expiry date</li> <li>• Do not connect to systems intended for plume evacuation, AGSS or vacuum</li> <li>• Do not try to disassemble the double mask</li> <li>• Do not use: Hypochlorite, phenol (&gt;5%), formalin, ketone, chlorinated hydrocarbons, aromatic hydrocarbons, inorganic acids, or denatured ethanol with concentrations higher than 70%</li> <li>• Sterilization with autoclave has not been validated</li> <li>• MRI compatibility has not been tested</li> </ul>	
<b>Limitations on reuse</b>	The double mask, coupling house and demand coupling house shall be replaced no later than 5 years from the manufacturing date or in the event of material damage. Check the manufacturing date according to figure 1, page 3.	

<b>Initial treatment at the point of use</b>	Not applicable
<b>Preparation before cleaning</b>	Not applicable
<b>Manual cleaning</b>	Not applicable
<b>Automated cleaning</b>	Use a certified medical washer-disinfector designed for washing and disinfecting medical tubes, anaesthesia equipment and respiratory accessories. <ol style="list-style-type: none"> <li>1. Place the double mask and coupling houses face down in a basket on a shelf in the dishwasher.</li> <li><b>WARNING!</b> Do not try to disassemble the double mask.</li> <li>2. Choose a washing programme that fulfils at least the following requirements:             <ul style="list-style-type: none"> <li>≥ 2 min. for rinsing with cold water</li> <li>≥ 5 min. dishwashing between 55°C - 95°C (158°F-203°F), with an enzymatic, alkaline detergent (pH ≤ 11).</li> <li>≥ 2 min. rinsing with cold water</li> <li>≥ 2 min. final rinsing with cold water</li> </ul> </li> </ol> <b>WARNING!</b> Do not use: Hypochlorite, phenol (>5%), formalin, ketone, chlorinated hydrocarbons, aromatic hydrocarbons, inorganic acids, or denatured ethanol with concentrations higher than 70%
<b>Disinfection</b>	Select the disinfection programme in the washer-disinfector. Ensure that the following parameters are met: <ul style="list-style-type: none"> <li>- Thermal rinsing with deionized water at 93°C (200°F) for 5 min.</li> </ul>
<b>Drying</b>	<ul style="list-style-type: none"> <li>- If possible, use the washer-disinfector drying function. Otherwise, dry in a cabinet designed for medical tubes, anaesthesia equipment and respiratory devices at 50°C-70°C (122°F- 158°F) for ≥ 20 minutes.</li> <li>- Drying can be completed manually by using a sterile lint-free wipe.</li> </ul> If none of the above alternatives are available, dry at room air by placing the devices on a drying rack.
<b>Maintenance, inspection and testing</b>	<ul style="list-style-type: none"> <li>- When all items have been cleaned and dried, lubricate the O-rings in the coupling houses by applying a thin layer of Fomblin lubricant with a cotton swab to improve seal ability and increase product lifespan.</li> <li>- A visual inspection should be made prior to every use to ensure that the surface is smooth, and no cracks or fissures are present. The double mask and coupling houses may become brittle over time and break into small pieces that may enter the patient's airway. <b>WARNING!</b> Do not use if the devices are cracked or damaged.</li> <li>- An inspection of the expiry date is recommended with every use. The manufacturing date can be found as two clock stamps on the devices. See figure 1, page 3. <b>WARNING!</b> Do not use after expiry date.</li> </ul>
<b>Packaging</b>	No particular requirements
<b>Sterilization</b>	Not applicable
<b>Storage</b>	Store at room temperature away from direct sunlight.

The instructions provided above have been validated by Medicvent as being capable of preparing the double mask and coupling houses for reuse. It remains the responsibility of the processor to ensure that the reprocessing is performed using the adequate equipment and materials to achieve the desired result. This requires verification and/or validation and routine monitoring of the process.

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Tillverkare: Medicvent AB      Metoder: Rengöring, Desinficering, Torkning